



**MEDISTER**  
Autoclave Class B  
Mediclave Series

Your Partner for Sterilization  
and Infection Control



**MEDICLAVE - Pre and Post Vacuum Autoclaves**

# Flash steam sterilizers for rapid sterilization with vacuum and printer



**Mediclave 12B(Flash)**



**Mediclave 35B / 50 B**



**Mediclave 17B / 21B**

- » Fast sterilization
- » Fully Automatic
- » 8 Preset Sterilization Programs
- » Ster. Temp.: 105-134C
- » Ster. Time : 4-60 min.
- » Dry Time : 1-25 min.
- » Vacuum Times : 1-10 times
- » B&D / HELIX / VACUUM Testing
- » Cleaning Program
- » Printer & USB Output Available
- » Big Stainless Steel Water Tanks
- » High-quality Bacteria Filter

## TECHNICAL DATA :

Model	Unit	Mediclave 12B(Flash)	Mediclave 17B	Mediclave 21B	Mediclave 35 B	Mediclave 50 B
Capacity	Ltrs.	12	17	21	35	50
Inner dimensions	mm	Ø 230 X 300D	Ø 230 X 410D	Ø 230 X 500D	Ø 325 X 410D	Ø 325 X 600D
Outer dimensions (WxDxH)	mm	570 X 500 X 450	570 X 610 X 450	570 X 775 X 450	665 X 705 X 545	665 X 895 X 545
Loading type	-	Front Loading	Front Loading	Front Loading	Front Loading	Front Loading
Heater power	kW	2 KW	2.4 KW	2.4 KW	2.8 KW	3.2KW
Trays	nos	2	2	2	3	3
Special Cycles	-	Vacuum / Helix / BD Test	Vacuum / Helix / BD Test	Vacuum / Helix / BD Test	Vacuum / Helix / BD Test	Vacuum / Helix / BD Test
Operating Temp. / Pressure	°C / kPA	121°C /134°C	121°C /134°C	121°C /134°C	121°C /134°C	121°C /134°C
Duty Cycle	-	8 Cycles / Day	8 Cycles / Day	8 Cycles / Day	8 Cycles / Day	8 Cycles / Day
Min Installation Size	mm	600 X 470	650 X 470	790 X 470	725 X 570	910 X 570
Total Cycle Time	min	20 min - 60 min	30 min - 75 min	30 min - 75 min	30 min - 75 min	30 min - 75 min
Cycle Record	-	Thermal Printer	Thermal Printer	Thermal Printer	Thermal Printer	Thermal Printer
Packing Dimensions (WxDxH)	mm	620 X 550 X 490	620 X 660 X 490	620 X 825 X 490	700 X 750 X 600	710 X 950 X 600
Net / Gross Weight	kg	45	52	55	65	85

Note : USB is optional Feature



**Office :** G-119, Gr. Floor, Shagun Arcade Premises, Near Dindoshi Bus Depot, Film City Road, Malad (East), Mumbai 40097

Contact : 9323233345

### Our Sales / Service Network :

Maharashtra / New Delhi / West Bengal / Tamilnadu / Karnataka / Andhra Pradesh / Kerala / Madhya Pradesh / Gujarat / Rajasthan / Punjab / Uttar Pradesh / Orissa / Assam / Himachal / Uttarakhand / Goa / Bihar / Jharkhand / Nagpur / Chattisgarh.